

KIRSCH Disc-type pen-recording thermometer

Disc-type pen-recording thermometer for installation in KIRSCH refrigerators and freezers.

It has a battery-operated quartz mechanism (independent of mains supply). The exchangeable recording disc is sufficient for a 7-day period. The temperature is recorded on the disc with a blue felt pen. The temperature sensor of the disc-type pen-recording thermometer is located in a refrigerant flask inside the refrigerator or freezer, so that the temperature of the appliance's contents is measured rather than the air temperature. The front window of the disc-type pen-recording thermometer can be locked.

Measuring range: -35° to $+15^{\circ}$ °C or -60° to $+25^{\circ}$ °C



Ill.: Disc-type pen-recording thermometer for installation in refrigerators and freezers.

If you have any questions, please do not hesitate to contact us. We will be happy to assist you.

Philipp Kirsch GmbH

Manufacturer of refrigeration technology for the health sector

Im Lossenfeld 14

77731 Willstätt | Germany

Tel.: +49 (0) 781 9227-0

Fax: +49 (0) 781 9227-200

info@kirsch-medical.com